

Leeds Forensic Systems

Your Comparison Microscope Source

Forensic Comparison Microscope Service

Maintain your labs ASCLD/ISO accreditation with regularly scheduled **maintenance, calibration** and **recertification** from Leeds Forensic Service. Our professional, factory trained technicians, will provide you with fast and reliable service on your Leeds Comparison Microscopes. We offer a complete service program, which includes field service as well as factory repairs at Leeds Forensic Service facility.

Leeds Forensic Service provides the following checks and adjustments:

- Mechanical functions and adjustments.
- Parfocality of all magnification positions.
- Parcentricity of all magnification positions.
- All lubrication points.
- Collimation of optical system.
- Cleaning of all external optics.
- Cleaning of internal optics as required
- Certification of magnification match, LCD, measurement scale and imaging software (if applicable).
- Examine illumination system.



Leeds Forensic Service Rates

- Preventive maintenance, calibration and recertification of one Leeds Forensic Comparison Microscope - \$900.
- Travel and related expenses are additional and will be quoted upon request.
- Standard hourly rate for Leeds Forensic Comparison Microscope Service, not including parts, is \$150 per hour.
- Multiple unit and annual contract pricing available upon request.



- Most repairs can be done on-site and will be quoted for customer's approval before work is performed.
- Loaner microscopes and equipment may be available.
- All of our Service Technicians are factory trained and certified to work on Leeds Forensic Comparison Microscope systems.
- **Maintain ASCLD/ISO requirements with regularly scheduled service, calibration and recertification with NIST traceable standards.**

Please call us to schedule service or repairs:
Leeds Forensic Systems, Inc.
800 Boone Avenue North
Minneapolis, MN 55427
WWW.LEEDSFORENSICS.COM

(763) 546-8575
(800) 444-5333
Fax: (763) 546-4369
forensics@leedsmicro.com